

Technical data

Standards

All specifications in ATPD: 'Ambient temperature and pressure dry'

bellavista was designed in accordance with guidance documents EN60601-1-4. Parameters endangering the patient are redundantly monitored by hardware safety components.

| Standard | Title |
|----------------|---|
| EN 60601-1 | Medical electrical equipment |
| EN 60601-1-2 | Electromagnetic compatibility |
| EN 60601-1-4 | Programmable electrical medical systems |
| EN 60601-1-6 | Usability |
| EN 60601-2-12 | Lung ventilators for medical use |
| EN ISO 10651-2 | Home care ventilators for ventilator-dependent patients |
| EN ISO 10651-6 | Home-care ventilatory support devices |
| EN ISO 14971 | Application of risk management to medical devices |

Classification

| Classification | Specification |
|-----------------------------------|---------------|
| Medical device class | 2b |
| Electrical protection class | Class I |
| Protection class of applied parts | Type BF |

Ventilator data

| Parameter | Specification |
|------------------------------------|--|
| Dimensions w x h x d | 35 cm x 22 cm x 33 cm |
| Weight | 8.5 kg |
| Internal batteries Battery life | Lithium ions 14.4 V / 6450 mAh ≈ 2 hours with one battery ≈ 4 hours with opt. second battery |
| Alarm volume | 45...86 dB(A) adjustable |
| Operational life span | 5 years |

Connection data

| Parameter | Specification |
|----------------------------|--|
| Power supply | 100...240 VAC, 50...60 Hz (80...264 VAC max. tolerance) |
| Power consumption | Typically: 80 W / max: 200 W |
| DC | 24 VDC (20...29 VDC max. tolerance) Type: 3.5 A / max: 7 A |
| Low-pressure oxygen supply | 0...150 kPa (0...1.5 bar) 0...15 L/min |
| USB connectors | 3x USB 2.0 |
| Ethernet connector | 1x 100 MBit |
| bellavista Bus | CAN Bus |
| Nurse call | 1) Closed on Alarm 2) Common 3) Open on Alarm |



Pressure, flow, resistance

| Parameter | Specification |
|--------------------------------|---|
| Max. possible pressure | 60 mbar |
| Max. flow @ 0 mbar | ≥ 260 L/min |
| Max. flow @ 60 mbar | ≥ 130 L/min |
| System resistance at 60 L/min | < 6 mbar/L/s |
| Generation of ventilation flow | Blower (turbine) and proportional valve for flow regulation |


Ambient conditions

| Parameter | Specification |
|--|-----------------------------|
| In operation | |
| Temperature | +5...+45 °C |
| Air pressure, altitude above sea level | 600...1100 hPa ≅ 4000 a.s.l |
| Relative humidity | 10...95 %RH, non condensing |
| Storage and transportation | |
| Temperature | -25...65 °C |
| Air pressure | 500...1100 hPa ≅ 4000 a.s.l |
| Relative humidity | 10...80 %RH, non condensing |

Supported ventilation modes

| Abbreviation | Description |
|--------------|---|
| CPAP | Continuous Positive Airway Pressure (p. Error! Bookmark not defined.) |
| P-A/C | Pressure Assist-Control Ventilation (p. Error! Bookmark not defined.) |
| PC-SIMV | Pressure Controlled - Synchronized Intermittent Mandatory Ventilation (p. Error! Bookmark not defined.) |
| PSV | Pressure Support Ventilation, (p. Error! Bookmark not defined.) |
| V-A/C | Volume Assist-Control Ventilation (p. Error! Bookmark not defined.) |
| VC-SIMV | Volume Controlled - Synchronized Intermittent Mandatory Ventilation (p. Error! Bookmark not defined.) |
| Day / Night | Automatic day / night switching (p. Error! Bookmark not defined.) |
| Dual Vent | Automatic switching between two modes (p. Error! Bookmark not defined.) |
| Mask Fit | Fitting the mask (p. Error! Bookmark not defined.) |
| Wean Vent | Weaning from the ventilator (p. Error! Bookmark not defined.) |
| Backup | Automatic ventilation when an apnoea occurs (p. Error! Bookmark not defined.) |
| Manual | Manually triggering a breath |
| Hold | Hold end-inspiratory or end-expiratory ventilation (≤ 30 s) |

Adjustable ventilation parameters

| Ventilation parameter | Adults | Children |
|---|---|--------------------------|
| P _{insp} [mbar] | 5...55 (12) | 3...45 (8) |
| P _{support} [mbar] | 5...55 (15) | 3...45 (10) |
| CPAP [mbar] | 4...20 (5) | 4...20 (4) |
| PEEP [mbar] | 0...40 (0) | 0...30 (4) |
| Rate AZ/min | 1...50 (12) | 1...100(20) |
| T _{insp} s | 0.2...10 (1.7) | 0.1...10 (1.3) |
| Trigger | Off / pressure / flow *) | Off / pressure / flow *) |
| Flow Trigger [L/min] | 0.1...15 (4) | 0.1...15 (2) |
| Pressure Trigger [mbar] | 0.1...10 (2) | 0.1...10 (1) |
| Exh Sens (expiration switching) [% max. insp. flow] | 10...90 (25) | |
| Ramp slope pressure increase [% inspiratory time] | 15...75 (15) | |
| V _{insp} [ml] | 300...2500 (500) | 40...500 (100) |
| Plateau [% of T _{insp}] | 0...70 | |
| Insp. flow curve |  | |

*) with optional proximal flow sensor
(Values in parentheses: default values)

Curves and monitoring parameters

| Parameter | Range | Resolution | Accuracy |
|------------------------|------------------|------------|------------------------------|
| Pressure ^{a)} | 0...100 mbar | 0.1 | ±2 mbar |
| Flow ^{a)} | -300...300 L/min | 0.1 | ±1 L/min ; 10% ^{b)} |
| Volume ^{a)} | 0...2500 mL | 1 | ±10 mL ; 10% ^{b)} |

| Parameter | Range | Resolution | Accuracy |
|---|----------------|------------|--|
| P _{Peak} | 0...100 mbar | 0.1 | ±4% ^{c)} |
| P _{Mean} | 0...100 mbar | 0.1 | ±4% ^{c)} |
| P _{Plateau} | 0...100 mbar | 0.1 | ±4% ^{c)} |
| PEEP | -10...100 mbar | 0.1 | ±4% ^{c)} |
| Rate | 0...200 AZ/min | 1 | ±1 |
| I:E | 1:10...10:1 | 0.1 | 10% ^{c)} |
| T _{insp} | 0...100 s | 0.1 | 10% ^{c)} |
| T _{Exp} | 0...100 s | 0.1 | 10% ^{c)} |
| V _{tinsp} | 0...999 mL | 0.1 | ±15% ^{c)} |
| V _{tExp} | 0...999 mL | 0.1 | ±15% ^{c)} |
| MV _{Exp} | 0...99.9 L/min | 0.01 | ±15% ^{c)} |
| FiO ₂ | 18...100 Vol% | 1 | ±(2.5%O ₂ + 2.5%) ^{d)} |
| SpO ₂ ^{*)} | 0...100 % | 1 | p. 5 |
| Puls ^{*)} | 18...300 1/min | 1 | p. 5 |
| Cardio Pleth ^{a)*)} | 0...255 | 1 | - |
| CO ₂ _{insp} ^{*)} | 0...25 Vol% | 0.1 | p. 5 |
| CO ₂ _{Exp} ^{*)} | 0...25 Vol% | 0.1 | p. 5 |
| CO ₂ ^{a)} | 0...25 Vol% | 0.1 | p. 5 |
| Time A | 99 days | 1s | 1s |
| Time B | 99 days | 1s | 1s |

*) With optional external sensor (p. **Error! Bookmark not defined.**)

a) Curve display

b) Larger value: specified value ; % of reading

c) % of reading

d) Absolute value + % of reading

Trending

| Parameter | Specification |
|---------------------------|--|
| Duration of the recording | 2 Weeks, after which the oldest data is automatically deleted. |
| What is recorded? | All curves and monitoring parameters |

Alarm limits

| Parameter | Adults | Children | Priority, Notes |
|--|--|---|------------------|
| V _{t_{exp}} [mL] *) | Δ : 300...3500 ∇ : 0...2500 | Δ : 40...700 ∇ : 0...500 | High Medium |
| P _{Peak} [mbar] | Δ : 0...60 | | High |
| MinVol [L/min] *) | Δ : 0.2...30 ∇ : 0.1...20 | | High Medium |
| Rate [AZ/min] | Δ : 1...100 ∇ : 0...99 | Δ : 1...150 ∇ : 0...149 | Medium Medium |
| Apnoea time [s] | Δ : 4...100 | Δ : 2...100 | High |
| FiO ₂ [%] | Δ : 21...100 ∇ : 19...99 | | Medium Medium |
| SpO ₂ [%] *) | ∇ : 0...99 | | |
| CO ₂ _{insp} [%] *) | Δ : 0.1...25 | | |
| CO ₂ _{Exp} [%] *) | Δ : 0.1...25 ∇ : 0...24.9 | | |

Δ Upper limit, ∇ lower limit

*) *) with optional proximal flow sensor, SpO₂ sensor or CO₂ sensor

Danger Do not use extreme alarm settings, as this could prevent the alarm from going off in case of an emergency.

Tubing system + sterilising filter

| Specification | Value |
|---|----------------------|
| Diameter | 22 mm |
| Standard | EN 12342 |
| Length with all accessories | 1...2 m |
| Max. flow resistance with all accessories | ≤4 mbar at 60 L/min |
| Sterilising filter | ≥ 99.999% efficiency |
| Flow resistance | <2 mbar at 60 L/min |

Only use our manufacturer recommended tubing systems (p. 4)

Do not use antistatic tubing systems!

SpO₂ finger clip oximeter (optional)

| Specification | Value |
|------------------------------|--|
| Dimensions | 32 x 32 x 51 mm cord 3 m |
| Measuring range and accuracy | 0...100% SpO ₂ ± 2% SpO ₂ ± 3% SpO ₂ ¹⁾ |
| | 18...300 pulse ± 3 1/min ± 5 1/min ²⁾ |
| Ambient conditions | *) |
| Measuring principle | Infrared at 600 nm and 910 nm |

1) When in motion or with poor circulation

2) When in motion, in the range of 40...210 pulse

*) bellavista specifications apply

CO₂ breathing gas sensor (optional)

| Specification | Value |
|--|--|
| Dimensions | 39 x 38 x 34 mm Cord 2.5 m |
| Measuring range and accuracy under standard conditions | 0...15 vol%CO ₂ ±(0.2 vol% + 2% of reading) |
| Measuring range and accuracy under all conditions | ±(0.3 vol% + 4% of reading) |
| Response time | <1 s ≤ 90 ms rise time 10...90% of reading |
| Ambient condition operation | *) |
| Ambient condition storage and transportation | -20...50 °C 500...1200 hPa 5...100 %RH condensing (dry for 24 h prior to use) |
| Warm-up time | 10s |
| Measuring principle | Absorption spectrum at 4.2 μm and 4.5 μm wave length |
| LED display on sensor | Green: sensor working Green flashing: Zero point calibration in progress Red: sensor defective Red flashing: airway adapter installed incorrectly |


*) For operation, bellavista specifications apply

Accessories, consumables, spare parts

| Art No | Description | Picture |
|----------------|---|---|
| 301.100.000 | bellavista 1000 ventilator completely packaged, with accessories |  |
| 301.105.000 | bellavista trolley |  |
| 301.106.000 | 2-flexible tubing arm for bellavista trolley |  |
| 301.113.000 | SpO ₂ pulseoximetry (p. Error! Bookmark not defined.) |  |
| 301.114.000 | CO ₂ capnography (p. Error! Bookmark not defined.) |  |
| 300.160.000 *) | Airway adapter for CO ₂ sensor (Box of 25 pieces) |  |
| 300.941.000 *) | Tubing system A Patient tubing system with inner pressure tubing |  |
| 301.167.000 *) | Tubing system C Patient tubing system for single use with expiratory valve | |
| 301.094.000 *) | Tubing system D Patient tubing system for single use with expiratory valve and flow sensor | |

| Art No | Description | Picture |
|----------------------------|--|---|
| 301.168.000 *) | Sterilising filter | |
| pend *) | flow sensor (single), comes with tubing system D | |
| 300.756.000 | EasyLung test lung |  |
| 300.964.000 | Leak adapter for use with test lung and tubing system A | |
| 301.108.000 301.109.000 | Power cord EUR 10 A, 2.5 m power cord CH 10 A, 2.5 m |  |
| 301.115.000 | Connector cable for curse call | |
| 301.116.000 | Connector cable for bellavista Bus | |
| 301.110.000 | Connector cable for 24 V feed | |
| 301.165.000 *) | Pack of filter mats (patient air filter + ventilator blower) (p. Error! Bookmark not defined.) | |
| 300.769.000 *) | O ₂ sensor (p. Error! Bookmark not defined.) |  |
| 300.999.000 | Fuse T 6.3 AH, 250 V (S. Error! Bookmark not defined.) |  |
| 301.097.000 | User manual bellavista 1000 | |
| 300.784.000 | Additional battery | |
| 301.158.000 | Accessory bag | |
| ?? | Annual maintenance | |

*) Consumable



The label contains technical information about bellavista, specific declarations and serves ventilator identification.













| Symbol | Explanation |
|-------------------------------------|---|
| | On/off switch |
| | Patient connector inspiration |
| i | Direct Access Button Provides information about connector options on the screen. |
| p | Connector for patient-centred pressure measurement. |
| | Connectors for patient-centred flow measurement Left: patient-centred pressure measurement Right: prox. flow sensor second tubing |
| | Connector for pressure-controlled expiratory valve. |
| BV-Bus | Connector for bellavista Bus |
| | Connector for nurse call |
| | Connectors for USB-devices |
| | Connector for Ethernet network connection |
| O2 In 1,5 bar max. | Connector for low-pressure oxygen supply Max: 1.5 bar / 15 L/min |
| DC In 24V/6A | Connector for external DC supply 24 V / 6 A |
| AC In 100 – 240V, 50 – 60 Hz | Connector for power supply 100...240 V, 50...60 Hz, 150 VA |

Symbols on ventilator and packaging


| Symbol | Explanation |
|-----------------------|---|
| | Manufacturer's specifications and date of manufacture |
| Schutzklasse I | The ventilator meets requirements of electrical protection class I |
| | Caution, note accompanying documents. |
| | All parts meet the requirements of protection class BF against electrical shock. |
| | The ventilator must not be disposed of with household waste. |
| CE 0124 | The CE symbol confirms conformity with the essential requirements of Directive 93/42/EU |
| | Temperature range for storage and transportation -25...65 °C |
| HR 10 – 80% | Humidity range for storage and transportation 10...80 %RH, non condensing |
| | Air pressure for storage and transportation 500...1100 hPa \approx 4000 a.s.l |
| | Fragile, please handle with care |
| | Transport in the upright position |
| | Recyclable packaging |
| | Keep dry |
| | Keep away from heat |

CO₂ breathing gas sensor (optional)

| Symbol | Explanation |
|--------|-------------|
|--------|-------------|

| Symbol | Explanation |
|---|---|
|  | Caution, note accompanying documents. |
|  | Model number |
|  | Lot number |
|  | Expiration date |
|  | Do not reuse |
|  | UL recognition symbol |
|  | Pressure control |
|  | Temperature control |
|  | Moisture control |
|  | Heat transfer |
|  | The CE symbol confirms conformity with the essential requirements of Directive 93/42/EU |
|  | In the US, by prescription only |

SpO₂ finger clip oximeter (optional)

| Symbol | Explanation |
|---|---|
|  | The CE symbol confirms conformity with the essential requirements of Directive 93/42/EU |
| SN | Serial number |

Declaration of conformity

imtmedical

Declaration of Conformity

The company
imtmedical ag
9470 Buchs, Switzerland

declares on sole responsibility that the production and the distribution of the following medical products of class IIb:


| Product | Type | |
|------------|--------------|---------|
| Bellavista | 1000 / 1000e | CE 0124 |


were made according to the quality rule EN ISO 9001 (quality management systems) in conjunction with EN ISO 13485 (quality assurance systems for medical products) and the guideline 93/42/ EWG appendix II.

This system meets all appropriate requirements for medical products of the guideline 93/42/ EWG issued by the Council of the European Community on 14th June 1993.

| Certificate center | Dekra Certification Services GmbH, Germany | | |
|-----------------------------------|--|---------------|-------------------------|
| Norm/ guideline | Date of re-certification | Valid to: | certificate registry no |
| Annex II, Section 3 von 93/42/EWG | March 10, 2006 | March 9, 2011 | 50747-16-01 |
| ISO 13485:2003 + AC:2007 | March 3, 2003 | March 9, 2011 | 50747-10-00 |
| ISO 9001:2008 | February 18, 2002 | March 9, 2011 | 50747-56-00 |

The conformity of the products is validated according to guideline 93/42/EWG, appendix II by the application of the sign

 *imtmedical ag*


Harri Friberg, CEO

Buchs, April 28, 2009

imtmedical

imtmedical ag Gewerbestrasse 8 9470 Buchs Switzerland
T +41 81 750 66 99 www.imtmedical.com www.bellavista-ventilator.com

All ® and ™ labeled brands and products are property of the company imtmedical and can not be used without written consent. Subject to changes.

301.182.000 | 11.09